

# CONTROL GROUP

## 1 Use case 1

ID: CPM:US1.001

Title: View vitals measurements

Actors: Actors: Doctor, nurse, social and health workers (only for home records), Caregivers (only for home records), Patient (only for home records)

Precondition:

Main scenario:

1. The user selects the 'View vitals measurements' tab for the given Patient.
2. The system searches for the measurements and displays the result list.
3. The actor browses the results.

Alternative flow no. 01:

1. If at step 3 of the main scenario, the actor selects the sub-option 'Add a measurement,' include the use case CPM:US1.002– Add a measurement.

Alternative flow no. 02:

1. If at step 3 of the main scenario, the actor selects the 'display trend,' include the use case CPM:US1.003– Display the graph of the measurements.

Exceptions:

### 1.1 Control group improved version

ID: CPM:US1.001

Title: View Vitals Measurements

Actors: **Doctor** can view the patient's recorded vitals measurements to monitor or change the treatment. **Caregivers**, **nurse**, and social-health **workers** can view the patient's recorded vitals measurements to monitor the Patient's health conditions and inform or alert the Patient's Doctor. **Patient** can view his (or her) vitals measurement to verify his (or her) health conditions.

**Preconditions:** At least one set of vitals measurements is stored for the given patient. **Vitals can be measured more than once per day.**

**Postconditions:**

- The user has successfully viewed the patient's vitals measurements.

Main scenario:

1. The user selects the "View Vitals Measurements" tab for the given patient.
2. The system searches for the measurements and displays them **in a table format, including the date and time of each measurement and the corresponding vital sign values.**
3. The **user views** the results.

Alternative flow no. 1:

1. If at step 3 of the main scenario, the actor selects the sub-option 'Add a measurement,' include the use case CPM:US1.002– Add a measurement.

# CONTROL GROUP

Alternative Flow no. 2:

1. If at step 3 of the main scenario, the actor selects the ‘display trend,’ include the use case CPM:US1.003– Display the graph of the measurements.

## Exceptions:

- If the system cannot find any measurements for the given patient, an error message is displayed to the user.
- If there is an error in displaying the measurements, an error message is displayed to the user.

## 2 Use Case 2

ID: CPM:US1.002

Title: Add a measurement

Actors: Doctor, nurse, patient, caregiver, social health operator

Precondition:

Main scenario:

1. The user clicks "New measurement".
2. The system displays a form in which it is possible to select the desired measuring instrument.
3. The user selects the desired measuring instrument.
4. The system searches and displays the vital parameters associated with that measuring instrument.
5. The user fills in the fields.
6. The user presses the “Save” button.
7. The system saves the changes for the new measurement to the list of measurements for the selected vital parameter.
8. The system displays the message "The measurement has been saved correctly".
9. The system displays the refreshed list of measurements with the new measurement.

Alternative flow no. 01:

1. At step 4 of the main scenario, the user presses the button “Cancel”.
2. The system closes the form.

Exceptions:

### 2.1 Control group improved version

ID: CPM:US1.002

Title: Add a measurement

Actors: **Doctor** can request the new measurement for treatment purpose. **Nurse**, **caregiver**, and **Social**-health operator can take the measurement for a given **Patient** and report to the Patient’s doctor. Patient can add the new measurement.

## Preconditions:

- The user is logged in.
- The user has permission to add a new measurement.

# CONTROL GROUP

## Postconditions:

- The system saves the new measurement to the measurement table.
- The user can list the new measurement for the given vitals.

## Main scenario:

1. The user clicks "New measurement".
2. The system displays a form with a list of measuring instruments.
3. The user selects the desired measuring instrument from the list.
4. The system searches and displays the vital parameters associated with the selected measuring instrument.
5. The user fills in the measurement fields relating to the vital parameters associated with the selected instrument (for example, for blood pressure, systolic and diastolic expressed in mmHg, for weight, the weight values expressed in kilograms, etc.).
6. The user presses the "Save" button.
7. The system validates the input data.
8. The system saves the new measurement to the list of measurements for the selected vital parameter.
9. The system displays the message "The measurement has been saved correctly."
10. The system refreshes the list of measurements to include the new measurement.

## Alternative flow no. 01:

3. At step 3 of the main scenario, the user decides not to add a new measurement.
4. The user clicks the "Cancel" button.
5. The system returns to the home page.

## Exceptions:

- 8a. The system displays the error message "The data is invalid."
- 10a. The system displays the error message "A system error occurred."

## 3 Use case 3

ID: CPM:US1.003

Title: Display the graph of the measurements for a vital parameter

Actors: Doctor, nurse, Social health operator, caregiver, patient

Preconditions:

Main scenario:

1. The user selects "Display measurement graph for vital" for the given parameter.
2. The system searches the measurements for the vital parameter and displays the related graph, having the measures on the ordinates, and the date-time on the abscissas.

# CONTROL GROUP

## 3.1 Control group improved

ID: CPM:US1.003

Title: Display Measurement Graph for a Vital Parameter

Actors: Doctor, nurse, social health operator, caregiver, patient

### Preconditions:

- The user is logged in and has the permissions to view the measurement graph.
- There is at least a measurement for the given vital parameter for the given Patient.

### Postconditions:

- The system displays the graph of the measurements for the given vital parameter.

Main Scenario:

1. The user selects "Display measurement graph for vital" for the given vital parameter.
2. The system searches for the measurements of the given vital parameter.
3. The system displays a graph with the vital parameter measurements on the y-axis and the date-time on the x-axis.
4. The user views the graph.

### Alternative Flow no. 1:

- At step 4, the user closes the Form.
- The system returns to the previous screen.

### Error Handling:

- If there is an error in generating the graph, the system displays an error message to the user and returns to the previous screen.
- If there are no measurements available for the selected vital parameter, the system displays an error message to the user.

## 4 Use case 4

ID: CPM:US1.004

Title: Add an automatic measurement

Actors: Patient, Sensor Device, Mobile App

Preconditions:

Main scenario:

1. The user uses an instrument from the kit to measure a vital sign.
2. The Sensor Device sends the raw data to the Mobile.
3. The Mobile App requests the System to save the measurement data.
4. The system saves the measurement.
5. The system sends the Mobile App the envelope data of the operation.

# CONTROL GROUP

6. The mobile App displays a message to the user to inform him (or her) that the data have been saved.

## 4.1 Control group improved version

ID: CPM:US1.004

Title: Add an automatic measurement

Actors: Patient, Sensor Device, Mobile App

### Preconditions:

- The patient is logged in and has the permissions to access the automatic measurement feature.
- The patient has a certified sensor device.
- The Mobile app is connected to Internet.

### Postconditions:

- The measurement data are saved for the given patient.
- The Mobile App displays an error message if the system was unable to insert the new data.

Main scenario:

1. The user selects a vital sign to measure using an instrument from the kit.
2. The Sensor Device sends the raw data to the Mobile App.
3. The Mobile App requests the system to save the measurement data (i.e. Vital sign type, date and time of measurement and values of the measurement)
4. The system validates the measurement data and saves them.
5. The system sends the Mobile App an envelope data of the operation, stating that the data were successful saved.
6. The Mobile App displays a message to the user to inform him (or her) that the data has been saved.

### Alternative flow no. 01:

- 4a. The user closes the app.
- 5a. The system does not save any data.

### Exception flow no. 01:

5a. The system is unable to save the measurement, and sends an envelope indicating the error to the Mobile App.

6a. The Mobile App displays an error message to the user stating that the measurement was not saved.

### Exception flow no. 2:

# CONTROL GROUP

2a. The device is not able to send data.

3a. The Mobile App displays an error message to the user stating that an error occurred with the device.

## Exception flow no. 3:

3a. The Mobile App is unable to connect to the system.

4a. The Mobile App displays an error message to the user stating that an error occurred.

## 5 Use case 5

ID: CPM:US2.001

Title: List all alert parameters

Actors: Doctor

Preconditions:

- The user selects the "Thresholds" in a patient's medical record

Main scenario:

1. The user selects the "List all vitals" button.
2. The system searches and displays the list of vitals.
3. The user selects the desired vital sign.
4. The system searches and displays the list of alerts related to the given vital parameter. The list is sortable by "Name".

### 5.1 Control group improved version

ID: CPM:US2.001

Title: List all alert parameters

Actors: Doctor

## Preconditions:

- The user must be logged in.
- The user must have permission to list alert parameters.
- The user has selected the "Thresholds" option.

## Postconditions:

- The user views the list of alerts for the given vital parameter.

Basic Flow:

1. The user selects the "List all vitals" button.
2. The system searches and displays the list of vitals.

# CONTROL GROUP

3. The user selects the desired vital sign.
4. The system searches and displays the list of alerts related to the selected vital parameter. The list is sorted by "Name" and contains the threshold value, the alert message, and the severity level.
5. The user can browse the results.

## Alternative Flow 1:

- 5a. The user clicks the “Close” button.
- 6a. The system closes the form and returns to the previous screen.

## Error Handling:

- 4a. The system displays the error message “A system error occurred.”
- 4b. The system displays the error message “No alerts are available for the selected vital sign.”